

Rhino Rescue


COMPRESSED GAUZE & ELASTIC CREPE BANDAGE



01 CUT the clothes of the injured part with scissors.



02 PREPARE compressed gauze to fill the wound.



03 EXTEND the casualty's arm at a 90-degree angle by placing it on your shoulder.



04 PACK the wound tightly with compressed gauze until the wound cavity is filled.



05 ENSURE the gauze extends 1-2 in above the skin.



06 HOLD pressure for a minimum of 3 min to ensure that the wound does not bleed.

Rhino Rescue

COMPRESSED GAUZE & ELASTIC CREPE BANDAGE



07 WRAP the elastic bandage around injured shoulder twice ensuring the gauze underneath is completely covered.



08 WRAP elastic bandage across, back and under the opposite armpit, anchoring around the opposite shoulder in a "figure 8" pattern.



09 After bandaging, SECURE with tie the tails of elastic bandage.



10 SECURE elastic bandage tails and knot using tape, prevent the knots from unraveling.

NOTE

The best position to treat the casualty is the seated position. If the casualty cannot be treated in the seated position, you will need to sit the casualty up as much as possible to apply the elastic bandage.